NEW

BINAPS®

Nasal Airflow Pressure and Snore Transducer

Monitor both Airflow and Snore via Salter Labs[®] Sleep Diagnostic Cannulas

The Salter Labs® BiNAPS® is a diagnostic pressure transducer device for use during sleep disorder studies. It will receive airflow/ snore pressure signals via Salter Labs® Sleep Diagnostic Oral/Nasal Cannulas, process through a special piezo type pressure transducer and deliver an accurate signal to your polysomnograph input box. Both airflow pressure and snore channels feature Hi – Lo selectors. They enable you to select the optimum signal for your recorder. The exclusive integral lead verification tester has a green light to indicate functional leads from the BiNAPS® to your input box. This is a great time saver to indicate appropriate signal transmission before beginning a study.

For Sleep Diagnostic Studies



Verify lead integrity before beginning the study. Don't miss another signal or have to repeat a test due to faulty or defective leads.

- Detect Airflow Pressure and Snore changes
- Hi/Lo switches to optimize airflow and snore signal levels into polygraph AC input box
- · No battery required
- Compact and lightweight
- Integral circuit tester/indicator for lead verification
- Self adhering hook and loop mounting patch included
- Lifetime warranty when used with Salter Labs[®]
 Sleep Diagnostic Cannulas

U.S. and Worldwide Patents Pending



BiNAPS® Nasal Airflow Pressure and Snore Transducer



Specifications

Sensor Technology	Piezo pressure transducer			
Output Types	Airflow and snore			
Size	5.75 cm	Width 3.48 cm 1.37"		
Weight	1.06 oz. / 30 grams			
Power Source	No battery required			
Output Connectors	1.5 mm touch proof plug (female)			
Airflow Nasal Cannula Input Connector				
Input Pressure Range	±25 cm H2O			
Signal Output Levels	1000 μV peak to peak			
Operating Temperature	5° C (40°F) - 40° C (104° F)			
Operating/Storage Humidity	15 – 95% Non-condensing			
Test Lead Recepticle	1.5 mm plug (male)			

Ordering	Information

Units Part #

BiNAPS® Nasal Airflow Pressure/Snore Transducer

5500



Fax: 661.854.3850 • Toll Free Fax: 1.800.628.4690

Email: salterlabs@us.salterlabs.com • www.salterlabs.com