Inductive bands



Sensors that make sense



SleepSense Inductive Effort Bands and Interface Cables

Made of a high-quality, soft and flexible material, these bands provide great comfort to the patient.

Velcro tabs allow the bands to fit a wide range of sizes.

The generous width of each band minimizes the sliding and shifting that occurs during the night.

System specific interface cables allow direct replacement of Piezo sensors, or easy connection with any system that accepts inductive sensors.


Sleepsense Inductive Effort Bands share these advantages:

- Expected life of more than 40 uses low cost per use
- Latex free
- 5 color-coded sizes, each fitting a wide range of patients
- Compatible with most inductive monitors
- Washable, steriliziable, autoclavable
- Very light and comfortable for the patient
- Body Position Sensors can be easily applied to the band

SLP is certified to IS09001:2000 and IS013485:2003

In the USA, FDA restricts these devices for sale by or on the order of a physician





Our full range of sensors can be found on our website at www.sleepsense.com

Inductive Bands

SleepSense Inductive Bands are available in five different sizes:

9002-L40	Inductive Effort Band - Small (40 cm / 15 Inch)	
9002-L60	Inductive Effort Band - Medium (60 cm / 23 Inch)	
9002-L90	Inductive Effort Band - Large (90 cm / 35 Inch)	
9002-L120	Inductive Effort Band - X-Large (120 cm / 47 Inch)	
9002-L133	Inductive Effort Band - XX-Large (133 cm / 52 Inch)	

Piezo Replacement Cable

This cable allows you to connect inductive bands to any system that accepts standard safety connectors. A plug and play replacement for any Piezo Sensor. No adjustments or switches required. Internal battery good for over 700 nights.

9101	Inductive Interface Cable, Abdomen/ DIN safety connectors	
9102	Inductive Interface Cable, Thorax/ DIN safety connectors	
9103	Inductive Interface Cable, Abdomen/ Key connector	



Interface Cables

For specific systems

9101 / 9102 A3	Inductive Interface Cables (Alice 3 [®])	
9101 / 9102 A4	Inductive Interface Cables (Alice 4 [®])	
9101 / 9102 A5	Inductive Interface Cables (Alice 5 [®])	
9101 / 9102 PT	Inductive Interface Cables (Sandman™ Pocket™)	
PS 008	Inductive Interface Cable (Medcare / Embla systems)	



300 · 004